

Rapid loop-mediated isothermal amplification assays for grapevine yellows phytoplasmas on crude leaf-vein homogenate has the same performance as qPCR

Polona Kogovšek · Nataša Mehle · Anja Pugelj · Tjaša Jakomin ·
Hans-Josef Schroers · Maja Ravnikar · Marina Dermastia

Accepted: 19 October 2016 / Published online: 27 October 2016
© The Author(s) 2016. This article is published with open access at Springerlink.com

Abstract A fluorescence-based real-time loop-mediated isothermal amplification (LAMP) assay for ‘*Candidatus Phytoplasma solani*’ (Bois noir phytoplasma; BNp) detection was developed and optimised for rapid laboratory and on-site BNp detection. This assay is highly specific, rapid and as sensitive as qPCR. It was validated according to European and Mediterranean Plant Protection Organisation recommendations. In addition, 286 grapevine leaf samples from the 2015 growing season were tested with this new real-time LAMP assay and an assay previously developed for detection of Flavescence dorée phytoplasma (FDp). These LAMP assays for detection of both BNp and FDp used without any DNA extraction step, which is a required step for qPCR analysis, were comparably effective to qPCR, and positive results were obtained in less than 35 min.

Keywords Real-time LAMP · Grapevine yellows phytoplasma · Validation

Introduction

In grapevine (*Vitis vinifera* L.), phytoplasmas are associated with grapevine yellows (GYs) diseases, that occur in the majority of grapevine growing countries worldwide (Constable et al. 2003). ‘*Candidatus Phytoplasma solani*’ is from the stolbur group/ 16SrXII-A of phytoplasmas (Quaglino et al. 2013), and is associated with Bois noir disease (i.e., Bois noir phytoplasma; BNp). In Europe, Bois noir disease is widespread. On the other hand, Flavescence dorée is associated with phytoplasmas belonging to the 16SrV group (i.e., Flavescence dorée phytoplasma; FDp) and this is the most severe of the GYs diseases. As such, FDp is listed in the EU2000/29 Council Directive on Harmful Organisms, and in the European and Mediterranean Plant Protection Organisation (EPPO) A2 quarantine list of pests. Typical symptoms of all GYs include leaf curling and discoloration of leaf veins and lamina, inter-vein yellowing or reddening (according to the variety), uneven or total lack of lignification of canes, flower abortion, and berry withering. Symptoms caused by BNp and FDp are not distinguishable by visual inspection and as the distribution of phytoplasma is uneven within a host that has a very low titre (Prezelj et al. 2013), only specific molecular approaches are suitable for the accurate and reliable detection of BNp and FDp.

Electronic supplementary material The online version of this article (doi:10.1007/s10658-016-1070-z) contains supplementary material, which is available to authorized users.

P. Kogovšek (✉) · N. Mehle · A. Pugelj · T. Jakomin ·
M. Ravnikar · M. Dermastia
Department of Biotechnology and Systems Biology, National
Institute of Biology, Večna pot 111, 1000 Ljubljana, Slovenia
e-mail: polona.kogovsek@nib.si

H.-J. Schroers
Agricultural Institute of Slovenia, Hacquetova ulica 17,
1000 Ljubljana, Slovenia

M. Ravnikar
University of Nova gorica, Vipavska cesta 13, 5000 Nova Gorica,
Slovenia

In diagnostic laboratories, assays based on real-time PCR (qPCR) are widely used for the detection of BNp and FDp (Angelini et al. 2007; Bianco et al. 2004; Galetto et al. 2005; Hren et al. 2007; Margaria et al. 2009; Pelletier et al. 2009). Furthermore, a protocol for the detection of FDp based on loop-mediated isothermal amplification (LAMP) assay was developed recently (Notomi et al. 2000; Kogovšek et al. 2015). LAMP is a highly specific and rapid technique, and it also circumvents the sensitivity of PCR and qPCR to inhibitors in plant extracts (Francois et al. 2011); furthermore, its isothermal nature provides the potential for it to be deployed in the field (Kogovšek et al. 2015; Tomlinson et al. 2010a). LAMP has shown a comparable or better performance to other detection methods and a wide applicability for the detection of plant pathogenic bacteria (Lenarčič et al. 2014), viroids (Lenarčič et al. 2013), fungi (Tomlinson et al. 2010b) and phytoplasmas (Bekele et al. 2011; Dickinson 2015; Hodgetts et al. 2011; Kogovšek et al. 2015; Tomlinson et al. 2010a).

In the present study, seven fluorescence-based real-time LAMP assays were designed for the specific detection of BNp and the best performing assay was verified. Additionally, a protocol was developed allowing rapid laboratory testing of phytoplasmas in grapevine crude leaf-vein homogenates. The performance of the LAMP assay was compared with the currently used qPCR protocols. This protocol includes the newly developed LAMP assay for BNp and the published LAMP assay for FDp (Kogovšek et al. 2015). These LAMP assays were shown here to be effective in the detection of both BNp and FDp in grapevine without any DNA extraction step, and with comparable sensitivity to that of qPCR. These LAMP protocols were also evaluated for on-site applications.

Materials and methods

Plant material, phytoplasma isolates and other pathogens

The grapevine samples used for testing and validation of the LAMP assay were collected from Slovenian vineyards during the 2013 and 2014 growing seasons (Table 1). In addition, for the comparative study of LAMP and qPCR, 286 grapevine samples were collected in 2015 during an official survey for BNp and FDp. Excised leaf veins and berries, berry pedicels and

tendrils were picked from symptomatic plants and stored at $-20\text{ }^{\circ}\text{C}$ for further analysis. The status of the plants (i.e., BNp/FDp infected, non-infected) was confirmed by qPCR (Hren et al. 2007) following the procedure described by Mehle et al. (2013a, b). The healthy and BNp-infected grapevine samples included in the present study were of different cultivars (Supplementary Table S1 and S5). The phytoplasma strains from different 16Sr groups that were included in this study originated from the collections of the National Institute of Biology (NIB, Ljubljana, Slovenia), University of Bologna (Italy), University of Nottingham (UK), French National Institute for Agricultural Research (INRA, France), and University of Udine (Italy) (Table 1). In addition, the tests included: (i) plant pathogenic bacteria from the NIB collections (including uncharacterised grapevine bacterial isolates), the National Collection of Plant Pathogenic Bacteria (York, UK) and the collection of Instituto Valenciano de Investigaciones Agrarias (IVIA, Spain), and (ii) several taxa of fungi frequently occurring in vineyards from the collection of the Agricultural Institute of Slovenia (KIS, Ljubljana, Slovenia) (Table 1).

LAMP primer design and reactions

The sequences for the rRNA genes (i.e., 16S, 23S, ITS region) and *secA*, *secY*, *stamp* and *tuf* were retrieved from the NCBI and aligned using the VectorNTI software (InforMax). The regions of homology specific for the stolbur phytoplasma strain (16SrXII-A) were determined. The LAMP assays were designed for all seven regions, using the LAMP Designer software (Premier Biosoft), and primers were synthesised by Integrated DNA Technologies (Germany).

All LAMP reactions were performed in a 25- μL reaction volume that contained 5 μL sample DNA or homogenate, 2 \times Isothermal Master Mix (Optigene), 0.2 μM F3 and B3 primers, 2 μM FIP and BIP primers, and 1 μM F-loop and B-loop primers. The LAMP reactions were performed in 8-well strips, or 96-well plate in GenieII (Optigene) or Roche LC480 instruments, respectively. In the ‘multi’ LAMP (mLAMP) assays where *secA* (targeting BNp; this study) and 23S rRNA (targeting FDp) (Kogovšek et al. 2015) were combined, the same primer and reagent concentrations were used as above. For the LAMP product annealing temperature determination (i.e., T_m), the fluorescence was detected during the cooling of the samples from 98 $^{\circ}\text{C}$ to 80 $^{\circ}\text{C}$

Table 1 Details of the samples used in the present study. The specificity of the *secA* LAMP assay was evaluated by testing different phytoplasma groups and subgroups isolated from different hosts, and bacteria and fungi isolates, and healthy grapevine samples. The sources of the phytoplasma-infected plant materials are given as the locations of the institutions of the individuals mentioned in the Acknowledgements. The phytoplasmas and bacteria from the field sampling were collected from 2006 to 2014, and were stored in the collection at the NIB

Host	Sample tissue type	Organism (acronym)	16Sr group (phytoplasma)	Material Source	Type	No. sample	LAMP	
							Tpos (min ± SD)	Tm (°C)
Phytoplasma								
<i>C. roseus</i>	Leaf	Brijnal little leaf (BLL)	16SH	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Rehmannia glutinosa (RG)	16SH-B	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Aster yellows (AYs)	16SH-B	Udine	Plant	1	-	-
<i>C. roseus</i>	Leaf	Strawberry green petal (SGP)	16SH-C	Nottingham	DNA	2	-	-
<i>C. roseus</i>	Leaf	Ribes in vinica (RIV)	16SH-C	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Vinca coconut phyllody (VCP)	16SHI	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Faba bean phyllody (FBP)	16SHI-C	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	<i>Crotalaria salhiana</i> phyllody (FBPSA)	16SHI-C	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Sweet potato little leaf (SPLL)	16SHI-D	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Western X-disease (WX)	16SHII-A	Bologna, INRA	DNA	7	-	-
<i>C. roseus</i>	Leaf	Lethal yellowing disease (LYD)	16SHV	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Elm yellows (EY-1)	16SRV-A	Udine	Plant	1	-	-
<i>C. roseus</i>	Leaf	FD70	16SRV	INRA	DNA	1	-	-
<i>V. vinifera</i>	Leaf	FD-D	16SRV	Field	Plant	4	-	-
<i>V. vinifera</i>	Leaf	FD-C	16SRV	Field	Plant	1	-	-
<i>V. vinifera</i>	Leaf	FD-C	16SRV	Field	Plant	1	-	-
<i>Clematis vitalba</i>	Leaf	FD-C	16SRV	Field	Plant	1	-	-
<i>C. roseus</i>	Leaf	Tomato big bud (TBB)	16SRVI	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Apple proliferation (API5)	16SPX-A	Udine	Plant	1	-	-
<i>C. roseus</i>	Leaf	European stonefruit yellows (ESFY)	16SPX-B	Udine	Plant	1	-	-
<i>Pyrus</i>	Leaf	Pear decline (PD)	16SPX-C	Field	Plant	1	-	-
<i>C. roseus</i>	Leaf	Napier grass stunt (NGS)	16SXI	Nottingham	DNA	1	-	-
<i>C. roseus</i>	Leaf	Stolbur (StolSE)	16SPXII-A	Udine	Plant	1	6.0	85.2
<i>C. roseus</i>	Leaf	BN 'tuf a'	16SPXII	AIT	DNA	1	8.9	85.2
<i>C. roseus</i>	Leaf	BN 'tuf b'	16SPXII	AIT	DNA	2	5.0 ± 0.0	85.1 ± 0.0
<i>V. vinifera</i>	Leaf	BN 'tuf a'	16SPXII	Field	Plant	6	6.2 ± 1.6	85.3 ± 0.1
<i>V. vinifera</i>	Leaf	BN 'tuf b'	16SPXII	Field	Plant	15	5.6 ± 0.9	85.3 ± 0.1
<i>V. vinifera</i>	Leaf	BN (not determined)	16SPXII	Field	Plant	36	12.0 ± 5.1	85.3 ± 0.4
Bacteria								
<i>V. vinifera</i>	Leaf	Unknown pathogenic bacteria, similar to <i>X. ampelinus</i>	Field	Plant	2	-	-	-
<i>V. vinifera</i>	Leaf	<i>Xylophilus ampelinus</i>	Field	Plant	3	-	-	-

Table 1 (continued)

Host	Sample tissue type	Organism (acronym)	16Sr group (phytoplasma)	Material		No. sample	LAMP		
				Source	Type		Tpos (min ± SD)	Tm (°C)	
<i>V. vinifera</i>		<i>Agrobacterium vitis</i>	KIS	(Av 13–2) ^a	1	-			
<i>V. vinifera</i>		<i>Agrobacterium vitis</i>	IVIA	(339–26) ^a	1	-			
Not determined		<i>Agrobacterium tumefaciens</i>	IVIA	(C58) ^a	1	-			
Not determined		<i>Agrobacterium rhizogenes</i>	IVIA	(K84) ^a	1	-			
<i>V. vinifera</i>	Leaf	Non-identified bacterial isolates	NIB ^b	Culture	22	-			
Fungi ^c									
<i>V. vinifera</i>	Leaf/berry	<i>Aspergillus jensenii</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Aureobasidium pullulans</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Bjerkandera</i> sp.	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Botrytis cinerea</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Cladosporium cf. cladosporioides</i>	Field	Plant	2	-			
<i>V. vinifera</i>	Leaf/berry	<i>Cryptococcus victoratae</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Trunk	<i>Diaporthe eres</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Didymella exitialis</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Trunk	<i>Fomitiporia mediterranea</i>	Field	Plant	2	-			
<i>V. vinifera</i>	Leaf/berry	<i>Fusarium cf. tricinatum</i>	Field	Plant	2	-			
<i>V. vinifera</i>	Leaf/berry	<i>Meyeromyces guilliermondii</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Penicillium brevicompactum</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Trunk, leaf/berry	<i>Phaeoaniella chlamydospora</i>	Field	Plant	4	-			
<i>V. vinifera</i>	Trunk	<i>Phanerochaete cf. chrysosporium</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Phlebia</i> sp.	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Schizophyllum</i> sp.	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Sporobolomyces roseus</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Thysanophora penicillitoides</i>	Field	Plant	1	-			
<i>V. vinifera</i>	Leaf/berry	<i>Trametes versicolor</i>	Field	Plant	1	-			
Healthy plant									
<i>V. vinifera</i>	Leaf	BNp-negative plants	Field	Plant	21	-			

BN(p), Bois noir (phytoplasma), Nottingham, University of Nottingham, UK; Bologna, University of Bologna, Italy; Udine, University of Udine, Italy; Field, from field sampling; INRA, French National Institute for Agricultural Research; AIT Austrian Institute of Technology; KIS, Agricultural Institute of Slovenia; IVIA, *Instituto Valenciano de Investigaciones Agrarias*, Spain; NIB, Slovenian National Institute of Biology; –, no amplification; *C. roseus*, *Catharanthus roseus*

^aBacterial collection; ^b Bacteria isolated from extracts of leaves from various grapevine cultivars grown on nutrient agar (Bacto Nutrient Agar, Difco) or yeast peptone glucose agar (yeast extract 7.0 g, proteose peptone 5.0 g, glucose 10.0 g, agar 15.0 g, distilled water to 1 L, pH 7.0); ^c Genomic DNA from pure cultured fungal isolates was extracted with the NucleoSpin Plant II (Macherey-Nagel) by using the cetyl-trimethylammonium-bromide-based extraction buffer supplied with this kit

(GenieII) or during the heating of the samples from 62 °C to 98 °C (Roche LC480). A sample was recognised as positive when the amplification curve was observed and the T_m was in the expected range in at least two replicates out of three. If only one out of three replicates was positive, the result was treated as inconclusive. All of the additional data regarding the Minimum Information for Publication of Quantitative Real-Time PCR Experiments guidelines (Bustin et al. 2010) for *secA* LAMP assay are given in Supplementary Table S2.

LAMP assay selection and optimisation

The performances of all of the seven LAMP assays were initially tested on a small set of BNp-positive, BNp-negative grapevine samples, and on phytoplasma isolates from other 16Sr groups, and bacteria. The LAMP assays were run at 62 °C and 65 °C. The LAMP assays that showed specific amplification and gave the shortest time of positivity (T_{pos}) at the selected temperature of amplification were taken forward for further testing (Supplementary Fig. S1).

Validation of the LAMP *secA* assay

The *secA* LAMP assay was validated according to the EPPO recommendations (EPPO 2014). The analytical sensitivity was evaluated by testing BNp-positive DNA that was diluted in DNA extracted from healthy grapevine midribs or in water. Three-fold dilutions were prepared and analysed with the LAMP and qPCR assays for the detection of BNp (Hren et al. 2007). This experiment was repeated four times, and each time the samples were analysed in triplicates, as a minimum, with each method. The analytical specificity was analysed by testing different BNp types (e.g., tuf a, tuf b), FDp, other phytoplasmas, bacteria, fungi and a healthy host (Table 1). Additionally, the sensitivities were compared to qPCR by testing a set of BNp-infected grapevine samples. The selectivity was evaluated by testing samples from various grapevine cultivars and tissues (Table 1). The repeatability was evaluated by analysing 3–5 replicates of DNA samples with various concentrations of BNp DNA. To assess reproducibility, the analyses were performed on up to 39 different days with freshly prepared reaction mix, and by up to three different operators and on two different devices (i.e., GenieII, Roche LC480), where one or more parameters, e.g., operator, device and/or day, were changed per repetition.

qPCR assays

All grapevine samples were tested for phytoplasmas (BNp, FDp) and the plant 18S rRNA gene (to confirm the DNA extraction success), with qPCR assays performed as previously described (Hren et al. 2007; Mehle et al. 2013b).

Leaf-vein crude homogenate testing using LAMP assays

The performance of *secA* LAMP assay was preliminarily evaluated on crude grapevine homogenates prepared with different homogenisation systems, i.e. different devices, lysing matrices, and buffers, thus assessing the applicability of different homogenisation approaches for direct testing (Supplementary Table S1). In a final comparative study we prepared samples with the FastPrep-24 (MP Biomedicals) protocol (40 s at 5 m/s), because of its comparably good performance with other protocols (Supplementary Table S1) and its routine use in a DNA extraction protocol (Mehle et al. 2013a). We additionally evaluated the analytical sensitivity of direct testing of leaf-vein homogenates prepared with Ultra Turrax Tube Drive device (UTTD, IKA). This device is favoured due to its cost and performance efficiencies and its on-site applicability for FDp detection (Kogovšek et al. 2015).

A comparative study was performed on 286 grapevine samples collected during an official survey in the 2015 growing season. One gram of grapevine plant material (i.e., leaf veins, tendrils, berries, berry pedicels) was added to 2 ml ELISA buffer (Kogovšek et al. 2015) in 15 mL test-tubes with garnet matrix and ceramic spheres, and homogenised with FastPrep-24 (MP Biomedicals) for 40 s at 5 m/s. The homogenates then underwent DNA extraction following the procedure described by Mehle et al. (2013a). The aliquots of the original homogenates were directly tested without prior DNA extraction. The homogenates and the extracted DNA were 10-fold diluted in water before their addition to the LAMP and qPCR assays, respectively. The crude grapevine leaf-vein homogenates were tested with the validated LAMP assay, targeting BNp, and 23S rRNA LAMP assay targeting FDp (Kogovšek et al. 2015) and the results compared to qPCR testing of extracted DNA. Additionally, the LAMP *COX* assay (Tomlinson et al. 2010b) was used for confirmation of the capability of amplifying plant DNA in the sample. The analytical

sensitivities of the on-site applicable leaf-vein crude homogenates testing were evaluated by testing five individual three-fold dilution series prepared from BNp-infected homogenates in healthy grapevine homogenate. The homogenates for the latter were prepared with an UTTD, and then directly tested with the LAMP assay. The results were compared to the standard BNp detection procedure in the laboratory that was based on qPCR analysis of the extracted DNA.

Results and discussion

In this study, seven LAMP assays were initially designed targeting diverse phytoplasma genes for the detection of BNp in grapevine samples. After the assay optimisation, their performances were compared. The best performing assay was validated according to the EPPO recommendations (EPPO 2014). In addition, the performance of the crude leaf-vein homogenate testing with the LAMP assay validated here for BNp and the earlier developed LAMP assay for FDp (Kogovšek et al. 2015) was compared to routine qPCR testing of the extracted DNA. The procedure of on-site testing of grapevine homogenates was developed and evaluated as well.

Primer design and evaluation of the LAMP assays

Sets of primers were designed according to the 16S rRNA, 23S rRNA, ITS, *secA*, *secY*, *stamp* and *tuf* genes (Table 2, Supplementary Table S3), where regions of the sequence specific for 16SrXII-A phytoplasmas were identified. In terms of the Tpos, specificity and sensitivity, the performance of each LAMP assay was evaluated by analysing a set of samples with different titres of

BNp, along with negative plant samples. Almost all of assays performed better at 62 °C than at 65 °C (data not shown). The *secY* assay did not show any amplification under the conditions tested, and the 16S rRNA and ITS assays showed non-specific amplification with unrelated phytoplasmas (data not shown); these were thus excluded from further analysis. The remaining LAMP assays, i.e. 23S rRNA, *secA*, *stamp* and *tuf*, were further evaluated and the *secA*, *stamp* and *tuf* assays showed higher specificities than the 23S rRNA assays (see Table 1 for the *secA* LAMP assay, data not shown for the other assays). Although the *stamp* assay was sufficiently specific, its Tpos was longer than that of the *secA* and *tuf* assays (data not shown), and it was thus excluded from further testing. The *secA* LAMP assay showed higher sensitivity than the *tuf* LAMP assay (data not shown), and was therefore selected for further validation here.

To optimise the procedure for GYs phytoplasma detection, an mLAMP assay for simultaneous detection of BNp and FDp was applied and evaluated, where the *secA* assay developed in this study and the earlier created 23S rRNA (Kogovšek et al. 2015) LAMP assays were combined. The mLAMP assay gave reliable results in terms of the amplification and Tm analysis (data not shown). However, the Tm analysis of the amplified product does not allow for differentiation between the two GYs phytoplasmas, as both have a similar Tm; therefore, the mLAMP assay was not evaluated further.

Sensitivity of the *secA* LAMP assay

The sensitivity of the *secA* LAMP assay and the qPCR assay for the detection of BNp (Hren et al. 2007) was compared. The testing of the dilution series showed that the *secA* LAMP assay detected as low as 9–27 BNp copies in a reaction, with all five replicates identified as positive (Table 3). Lower amounts of the BNp DNA in the reaction resulted in lower numbers of positive replicates. Nevertheless, results of the *secA* LAMP assay suggested that it had a similar sensitivity compared to the qPCR assay for BNp detection (Supplementary Table S4).

The *secA* LAMP assay was further evaluated by testing grapevine leaf-vein samples with different amounts of BNp DNA, as estimated using qPCR. The diagnostic sensitivity of the *secA* LAMP assay was 100 %, which was supported by the detection of BNp in all of the 57 samples tested (Table 1, Supplementary

Table 2 *secA* LAMP primer sets designed for detection of BNp

Primer	Primer sequence (5'–3')
<i>secA</i> -F3	TTCCACCAAATCTTTGAGCT
<i>secA</i> -B3	ACAATAGCTACCAATATGGCAG
<i>secA</i> -FIP	AATTAAGAGGACGTGCCGGTTCGTTT GTCTTCGCTGGAAA
<i>secA</i> -BIP	TCCTAAAACCGCCAAACCTCCGAAG AGGAAGTATATTCGCTT
<i>secA</i> -LoopF	GAGATCCTGGTTATTCTCGCTT
<i>secA</i> -LoopB	TTCAACAACGCCTTACCT

Table 3 Analytical sensitivities of the *secA* LAMP assays tested on extracted DNA in parallel with the qPCR. Three-fold dilution series of BNp-infected DNA were prepared in healthy grapevine DNA, and tested in at least triplicate. The mean Cq, Tpos and Tm (\pm standard deviation) are shown

Fold dilution	Estimated copy number	qPCR (Cq \pm SD)	LAMP	
			Tpos (min \pm SD)	Tm ($^{\circ}$ C \pm SD)
10	729–2187	31.8 \pm 0.1	6.1 \pm 0.2	85.6 \pm 0.1
30	243–729	33.2 \pm 0.2	8.1 \pm 0.8	85.4 \pm 0.1
90	81–243	35.4 \pm 1.1	9.8 \pm 1.2	85.3 \pm 0.1
270	27–81	36.4 \pm 0.3	11.9 \pm 2.5	85.3 \pm 0.0
810	9–27	37.3 \pm 0.4 ^a	13.3 \pm 2.4	85.5 \pm 0.1
2430	3–9	37.6 \pm 0.4 ^b	21.3 \pm 12.4 ^c	85.2 \pm 0.1
7290	1–3	38.0 ^d	17.1 ^d	85.6
21,870	0	neg	neg	/
65,610	0	neg	neg	/

^aFour out of five replicates were positive; ^bthree out of five replicates were positive; ^ctwo out of five replicates were positive; ^done out of five replicates was positive; neg, no amplification was detected; /, no Tm value available

Table S5). However, in five DNA samples with low BNp titres (Cq >34.5), BNp DNA was only confirmed after the undiluted sample DNA was re-analysed. These data showed again that the LAMP reaction is less sensitive to inhibition, and allow analysis of undiluted grapevine DNA samples (Kogovšek et al. 2015). Therefore, testing of undiluted DNA is recommended to avoid false-negative results. All of the positive results from LAMP reactions were observed in less than 35 min of amplification, with a mean Tpos of 9.3 ± 4.6 min. Rapidity of these reactions suggests the possibility of its use for rapid diagnostics.

Specificity and selectivity

The analytical specificity of the *secA* LAMP assay was first evaluated by *in-silico* analysis and showed high predicted specificity for the 16SrXII-A phytoplasmas, including BNp (data not shown).

The specificity of the *secA* LAMP assay was confirmed through the analysis of the samples listed in Table 1 and in Supplementary Table S5. All of the data obtained with the *secA* LAMP assay were in agreement with the qPCR data obtained according to Hren et al. (2007). The developed assay was not designed to discriminate tuf type-a and tuf type-b BNpd (Foissac et al.

2013), therefore no difference in amplification of these tuf-types was detected. No cross-reactivity of the *secA* LAMP assay was observed among the phytoplasma DNA from other 16Sr groups, including FDP from grapevine samples from several European winegrowing areas, and the pathogenic and non-pathogenic bacteria and fungi and the healthy host (Table 1).

The Tm for the specific amplicon ranged from 84 $^{\circ}$ C to 85 $^{\circ}$ C on the GenieII apparatus (data not shown), and from 85 $^{\circ}$ C to 86 $^{\circ}$ C on the Roche LC480 system (Table 1). This variation between the different instruments is expected and Tm should be experimentally determined for the individual instruments (Kogovšek et al. 2015). In the same five 10-fold-diluted DNA samples with low BNp titres (Cq >34.5) discussed above with respect to the sensitivity, up to 2 $^{\circ}$ C higher Tm were determined, while the undiluted DNA from the samples gave a Tm in the expected range (Supplementary Table S5).

Repeatability and reproducibility

Repeatability and reproducibility of the assays were evaluated by analysis of several replicates of DNA samples with various BNp DNA titres. When testing replicates of the same sample with high and medium titres of BNp DNA (i.e., qPCR Cq <33), the assay repeatability was 100 %, meaning that all of the replicates gave a positive result (Supplementary Table S4). At lower BNp concentrations, the detection of the BNp DNA using the *secA* LAMP assay varied, giving positive and negative results, and this can be attributed to stochastic effects in the target copy distribution in the replicates (Hren et al. 2007). These data were 100 % reproducible when tested with the different devices, and by different operators, on different days, and with different reaction mixes (Supplementary Table S6).

Comparative study of the analysis of leaf-vein crude homogenates with the LAMP assay and qPCR testing of DNA

To determine the performance of the validated LAMP assays for the BNp detection (this study) and FDP detection (Kogovšek et al. 2015), an extensive comparative study was set up in 2015. We analysed all grapevine samples received for GYs disease testing during the official national survey. A total of 286 samples were processed. Crude leaf-vein homogenates were tested

with the LAMP assays, and were further subjected to DNA extraction for qPCR testing (Hren et al. 2007).

BNp was confirmed by qPCR BNp in 208 out of 286 samples and FDp in 24 samples, with four samples being positive for both BNp and FDp (Supplementary Table S8). With the LAMP assays, the sanitary status of 284 out of the 286 samples tested (i.e. 99.3 %) were correctly determined, and only two samples (D211/15, D333/15) were misinterpreted (Table 4). In sample D333/15, the concentration of BNp was low (mean Cq, 32.4), and it was even lower in sample D211/15, where BNp was detected with qPCR only in an undiluted DNA sample; testing of the 10-fold dilution of this DNA resulted in only one positive replicate (Cq, 38.9) out of three tested. In five BNp-infected samples with Cq >32, and in one FDp-infected sample with Cq 31.4, the results obtained by the LAMP assays were inconclusive. Only one positive replicate out of three tested indicated low phytoplasma levels in the sample. In the case of the rapid in-laboratory screening of the crude homogenates with LAMP assays, DNA would be extracted from these samples and the DNA tested by qPCR.

The performance of both of the LAMP assays on crude homogenates was also determined in terms of Tpos and Tm (Supplementary Table S8). All of the LAMP reactions gave positive results in up to 35 min of amplification, with means of 17.3 min and 12.7 min for FDp and BNp, respectively. The Tm values of the

Table 5 Analytical sensitivity of the *secA* LAMP assay applied in the direct homogenate testing in parallel with the qPCR for the extracted DNA. A three-fold dilution series of BNp-infected plant material was prepared in healthy grapevine material and tested in triplicate, except for the 270-fold diluted sample, which was tested in duplicate with the LAMP assay. The mean Cq, Tpos and Tm (\pm standard deviation) are shown

Fold dilution	Estimated copy number	qPCR (Cq \pm SD)	LAMP	
			Tpos (min; \pm SD)	Tm ($^{\circ}$ C; \pm SD)
10	81–243	32.4 \pm 0.1	20.4 \pm 3.5	85.5 \pm 0.1
30	27–81	34.0 \pm 0.2	14.7 \pm 3.2	85.4 \pm 0.1
90	9–27	35.0 \pm 0.5	24.8 \pm 14.1	85.2 \pm 0.1
270	3–9	36.0 \pm 0.4	15.4 ^a	85.2 ^a
810	1–3	36.2 \pm 0.8	neg	/
2430	0	neg	neg	/
7290	0	neg	neg	/
21,870	0	neg	neg	/

^a One replicate out of two tested was positive; neg, no amplification was detected; /, no Tm value available

amplified products were in the expected range (23S rRNA assay Tm, from 84 $^{\circ}$ C to 86 $^{\circ}$ C). The variation in the Tm observed for the *secA* LAMP assay (84 $^{\circ}$ C to 87 $^{\circ}$ C) was due to the non-optimal reaction conditions in one of the LAMP runs (samples D267/15–D275/15), confirmed by the irregularities in the shapes of the amplification curves (data not shown). However, not

Table 4 Overview of the LAMP and qPCR comparative study results

qPCR results (Cq)	No. of samples	LAMP = qPCR	LAMP \approx qPCR ^a	LAMP \neq qPCR
BNp neg, FDp neg	50	50	0	0
BNp neg, FDp pos (Cq <31)	22	22	0	0
BNp neg, FDp pos (Cq >31)	2	2	0	0
BNp pos (Cq >31), FDp pos (Cq <31)	1	0	1 ^b	0
BNp pos (Cq <31), FDp pos (Cq <31)	1	1	0	0
BNp pos (Cq <31), FDp pos (Cq >31)	2	1	1 ^c	0
BNp pos (Cq <31), FDp neg	187	187	0	0
BNp pos (Cq >31), FDp neg	21	15	4	2
Total	286	278	6	2

LAMP qPCR Number of samples equally diagnosed, LAMP qPCR Number of samples with inconclusive results with LAMP, LAMP \neq qPCR Number of samples differentially diagnosed. Neg Negative, pos Positive

^a One parallel out of three tested (1/3) was positive with LAMP

^b LAMP BN (1/3), LAMP FD positive

^c LAMP BN pos, LAMP FD (1/3)

including the T_m values of these samples, the T_m of the *secA* LAMP assays were from 84 °C to 85 °C. In both the *secA* and 23S rRNA LAMP assays, the T_m of the reactions on the crude grapevine leaf-vein homogenates was slightly lower than that in the reactions with the extracted DNA (Kogovšek et al. 2015). The variation in the T_m of the amplified product is most probably due to matrix interference (i.e., high plant-cell components) influencing the melting of the amplified product, which was previously observed at the DNA level. An influence of the sample matrix on T_m was also observed in the testing of the homogenates prepared with Na-acetate buffer, where the T_m varied by up to 2 °C in individual replicates of the same sample (data not shown). However, the exact mechanism(s) of the matrix influence on the melting curve analysis and consequent T_m values is not yet known.

In addition, homogenisation with the on-site applicable UTDD was further evaluated in terms of the analytical sensitivity. The homogenisation procedure with the UTDD was shown to allow the detection of as low as three to nine copies of BNp DNA in a homogenate sample (Table 5). Lower amounts of BNp DNA in the reaction reduced the number of positive replicates. Nevertheless, based on all of the data from the on-site BNp testing approach, it was similar to the qPCR assay, where extracted DNA was tested (Supplementary Table S7).

Conclusions

The *secA* LAMP assay was developed and validated for the detection of BNp in extracted grapevine DNA. The results showed that the assay is fast, highly specific, and similarly sensitive as qPCR. A new assay was validated in accordance with EPPPO recommendations. The assay can be used together with the LAMP assay for FDP detection (Kogovšek et al. 2015) and is innovative in-laboratory and on-site rapid testing of grapevine crude homogenates of different tissues. LAMP assays for BNp and FDP gave positive reaction results in 6–35 min which is therefore sufficient amplification time frame for GY phytoplasma detection. Performance of LAMP assays of selected grapevine yellows phytoplasmas on crude leaf vein homogenates provided similar results as the more laborious qPCR assay on extracted DNA from the same samples.

Acknowledgments The authors thank Ana Kompan, Larisa Gregur and Špela Prijatelj Novak for excellent technical support. Samples were taken during the official survey for BNp and FDP in Slovenia by the Slovenian Phytosanitary Administration. The authors would also like to thank Phil Jones (Rothamsted Research, UK), Assunta Bertaccini (University of Bologna, Italy), Elisabeth Boudon-Padiou and Xavier Foissac (INRA, France), Ruggero Osler and Paolo Ermacora (University of Udine, Italy), Matthew Dickinson (University of Nottingham, UK), Günter Brader (AIT Austrian Institute of Technology GmbH, Austria), *Instituto Valenciano de Investigaciones Agrarias* (IVIA, Spain), Slovenian Phytosanitary Administration Inspectors, and Tanja Dreo and Manca Pirc (NIB, Slovenia) for providing bacterial isolates/ samples. Fungal taxa included in this study were isolated in the context of FP7 ERA-Net Project Core Organic II, [Vineman.org](http://www.vinemanager.org). Language revision was carried out by Christopher Berrie.

Compliance with ethical standards

Conflict of interest The authors declare that they have no conflict of interest.

Funding This study was financially supported by the Slovenian Research Agency (project no. L4–5525) and by Euphresco Grafdepi 2.

Open Access This article is distributed under the terms of the Creative Commons Attribution 4.0 International License (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium, provided you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made.

References

- Angelini, E., Luca Bianchi, G., Filippin, L., Morassutti, C., & Borgo, M. (2007). A new TaqMan method for the identification of phytoplasmas associated with grapevine yellows by real-time PCR assay. *Journal of Microbiological Methods*, 68(3), 613–622.
- Bekele, B., Hodgetts, J., Tomlinson, J., Boonham, N., Nikolić, P., Swarbrick, P., et al. (2011). Use of a real-time LAMP isothermal assay for detecting 16SrII and XII phytoplasmas in fruit and weeds of the Ethiopian Rift Valley. *Plant Pathology*, 60(2), 345–355.
- Bianco, P. A., Casati, P., & Marziliano, N. (2004). Detection of phytoplasmas associated with grapevine flavescence doree disease using real-time PCR. *Journal of Plant Pathology*, 86, 257–261.
- Bustin, S. A., Beaulieu, J.-F., Huggett, J., Jaggi, R., Kibenge, F. S. B., Olsvik, P. A., et al. (2010). MIQE précis: practical implementation of minimum standard guidelines for fluorescence-based quantitative real-time PCR experiments. *BMC Molecular Biology*, 11. doi:10.1186/1471-2199-11-74.

- Constable, F. E., Gibb, K. S., & Symons, R. H. (2003). Seasonal distribution of phytoplasmas in Australian grapevines. *Plant Pathology*, *52*, 267–276.
- Dickinson, M. (2015). Loop-mediated isothermal amplification (LAMP) for detection of phytoplasmas in the field. In C. Lacomme (Ed.), *Methods in molecular biology: plant pathology* (Vol. 1302, pp. 99–111). New York: Springer.
- EPPO (2014). Specific requirements for laboratories preparing accreditation for a plant pest diagnostic activity PM7/98 (2). *EPPO Bulletin*, *44*, 117–147.
- Foissac, X., Carle, P., Fabre, A., Salar, P., Danet, J.-L., & Stolbur-Euromed consortium, NN (2013). “Candidatus phytoplasma solani” genome project and genetic diversity in the Euro-Mediterranean basin. 3rd European Bois noir Workshop (pp. 11–13). Barcelona: Proceedings of the 3rd European Bois noir Workshop.
- Francois, P., Tangomo, M., Hibbs, J., Bonetti, E. J., Boehme, C. C., Notomi, T., et al. (2011). Robustness of a loop-mediated isothermal amplification reaction for diagnostic applications. *FEMS Immunology and Medical Microbiology*, *62*, 41–48.
- Galetto, L., Bosco, D., & Marzachi, C. (2005). Universal and group-specific real-time PCR diagnosis of flavescence dorée (16Sr-V), bois noir (16Sr-XII) and apple proliferation (16Sr-X) phytoplasmas from field-collected plant hosts and insect vectors. *Annals of Applied Biology*, *147*(2), 191–201.
- Hodgetts, J., Tomlinson, J., Boonham, N., Gonzalez-Martin, I., Nikolić, P., Swarbrick, P., et al. (2011). Development of rapid in-field loop-mediated isothermal amplification (LAMP) assays for phytoplasmas. *Bulletin of Insectology*, *64*, 41–42.
- Hren, M., Boben, J., Rotter, A., Kralj, P., Gruden, K., & Ravnikar, M. (2007). Real-time PCR detection systems for flavescence dorée and bois noir phytoplasmas in grapevine: comparison with conventional PCR detection and application in diagnostics. *Plant Pathology*, *56*, 785–796.
- Kogovšek, P., Hodgetts, J., Hall, J., Prezelj, N., Nikolić, P., Mehle, N., et al. (2015). LAMP assay and rapid sample preparation method for on-site detection of flavescence dorée phytoplasma in grapevine. *Plant Pathology*, *64*(2), 286–296.
- Lenarčič, R., Morisset, D., Mehle, N., & Ravnikar, M. (2013). Fast real-time detection of potato spindle tuber viroid by RT-LAMP. *Plant Pathology*, *62*(5), 1147–1156.
- Lenarčič, R., Morisset, D., Pirc, M., Llop, P., Ravnikar, M., & Dreo, T. (2014). Loop-mediated isothermal amplification of specific endoglucanase gene sequence for detection of the bacterial wilt pathogen *Ralstonia solanacearum*. *PLoS One*, *9*(4), e96027.
- Margaria, P., Turina, M., & Palmano, S. (2009). Detection of flavescence dorée and bois noir phytoplasmas, *Grapevine leafroll associated virus-1* and *-3* and *Grapevine virus A* from the same crude extract by reverse-transcription real-time Taqman assays. *Plant Pathology*, *58*(5), 838–845.
- Mehle, N., Nikolić, P., Rupar, M., Boben, J., Ravnikar, M., & Dermastia, M. (2013a). Automated DNA extraction for large numbers of plant samples. In M. Dickinson & J. Hodgetts (Eds.), *Phytoplasma methods and protocols* (pp. 139–145). London: Humana Press.
- Mehle, N., Prezelj, N., Hren, M., Boben, J., Gruden, K., Ravnikar, M., et al. (2013b). A real-time PCR detection system for the bois noir and flavescence dorée phytoplasmas and quantification of the target DNA. In M. Dickinson & J. Hodgetts (Eds.), *Phytoplasma methods and protocols* (pp. 253–268). London: Humana Press.
- Notomi, T., Okayama, H., Masubuchi, H., Yonekawa, T., Watanabe, K., Amino, N., et al. (2000). Loop-mediated isothermal amplification of DNA. *Nucleic Acids Research*, *28*, e63.
- Pelletier, C., Salar, P., Gillet, J., Cloquemin, G., Very, P., Foissac, X., et al. (2009). Triplex real-time PCR assay for sensitive and simultaneous detection of grapevine phytoplasmas of the 16SrV and 16SrXII-A groups with an endogenous analytical control. *Vitis*, *48*(2), 87–95.
- Prezelj, N., Nikolić, P., Gruden, K., Ravnikar, M., & Dermastia, M. (2013). Spatiotemporal distribution of flavescence dorée phytoplasma in grapevine. *Plant Pathology*, *62*(4), 760–766.
- Quaglino, F., Zhao, Y., Casati, P., Bulgari, D., Bianco, P. A., Wei, W., et al. (2013). “Candidatus phytoplasma solani”, a novel taxon associated with stolbur- and bois noir-related diseases of plants. *International Journal of Systematic and Evolutionary Microbiology*, *63*(8), 2879–2894.
- Tomlinson, J. A., Boonham, N., & Dickinson, M. (2010a). Development and evaluation of a one-hour DNA extraction and loop-mediated isothermal amplification assay for rapid detection of phytoplasmas. *Plant Pathology*, *59*(3), 465–471.
- Tomlinson, J. A., Dickinson, M. J., & Boonham, N. (2010b). Rapid detection of *Phytophthora ramorum* and *P. kernoviae* by two-minute DNA extraction followed by isothermal amplification and amplicon detection by generic lateral flow device. *Phytopathology*, *100*, 143–149.